LUVION®

Premium Babyproducts

Exact non-contact thermometer

ENGLISH - USER MANUAL

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Thank you for purchasing the Luvion Exact non-contact digital Infrared thermometer.

The Exact non-contact thermometer is a high quality product incorporating the latest technology and it has been tested in accordance with international standards. With its unique technology, the Luvion Exact non-contact thermometer will guarantee reliable readings. The thermometer will self-test every time it starts up to always ensure the accuracy of measurements.

The Exact non-contact thermometer can be used to measure the body temperature, and is suitable for people of all ages. The thermometer can be used to measure the temperature of milk as well.

Please read these instructions carefully before using this instrument and keep the instructions and the thermometer in a safe and secure place.

Contents

1. The advantages of the Luvion Exact non-contact thermometer	3
2. Necessary safety instructions	3
3. Product description	5
4. How to measure the forehead temperature	5
5. How to measure the milk temperature	5
6. Explanation of displays and operation instructions	6
7. Cleaning and storage	10
8. Technical specifications	10
9. After-sale service	12
10. Security type	12

1. The advantages of the Luvion Exact non-contact thermometer

Measurement within seconds

The innovative infrared technology allows measurement of forehead temperature in just a few seconds.

Accurate and reliable

Due to the unique probe assembly construction, the advanced infrared sensor and the complete calibration process, this thermometer will offer a very accurate and reliable temperature measurement.

Gentle and easy to use

- The non-contact thermometer is non-invasive. A measurement can be taken even while a child is sleeping.
- The unique design makes the non-contact thermometer easy to hold and use.
- The non-contact thermometer is less threatening to a child than a rectal thermometer.

Multiple reading recalls

Users can recall the last 20 readings when activating the memory mode, enabling efficient tracking of temperature history and variations.

Safe and hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.

Fever Indication feature


The Fever Indication feature helps you to interpret the meaning of the temperature faster thanks to the colour indicated display.

The screen displays green if there is no fever and a red alert for fever. Seven audible beeps indicate temperatures above 37.5°C.

2. Necessary safety instructions

- Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.
- Never immerse the Forehead Thermometer into water or other liquids (not waterproof). For cleaning and disinfecting please follow the instructions in the "Cleaning and Storage" section.

- Keep the instrument away from direct exposure to the sun and keep it in a dust-free, dry area at a temperature between 16° - 40°C.
- Do not use the thermometer if there are signs of damage on the measuring sensor or on the instrument itself. If damaged, do not attempt to repair the instrument!
- This non-contact thermometer consists of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock.
Do not twist the instrument or the measuring sensor.
- The thermometer is only suitable for the measurement of body temperature and the temperature of liquids like milk. The use of this thermometer is not intended as a substitute for consultation with your doctor. The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to these instructions, the accuracy of measurements will not be affected.
- The non-contact thermometer is supplied with 2 AAA batteries. Replace with 2 new AAA batteries when the flashing battery symbol appears on the LCD display. To change the batteries, slide open the battery cover as shown and remove batteries. When replacing the batteries be sure to align the batteries properly as indicated inside the battery compartment.
- Remove the battery from the instrument if it is not required for extended periods of time, in order to avoid damage to the thermometer resulting from a leaking battery.


⚠ WARNING:

- Use of this non-contact thermometer is not intended as a substitute for consultation with your physician or pediatrician.
- This thermometer is not waterproof! Please NEVER immerse into liquids!

3. Product description

- (1) Sensor
- (2) Body measuring button
- (3) Milk measuring button
- (4) LCD display
- (5) Battery cover

4. How to measure the body temperature

The Luvion Exact non-contact thermometer measures infrared energy emitted from the skin above the eyebrow area and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. Temperature readings obtained by scanning the skin above the eyebrow area will provide the greatest accuracy.

The Luvion Exact non-contact thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

5. How to measure the milk temperature

The non-contact thermometer measures infrared energy emitted from the feeding bottle, and it shall display the milk temperature by transition calculating and emission rate compensation. Please use the button on the side of the thermometer (3) for startup when measuring milk temperature.

To avoid an inaccurate measurement, note the following:

- Patients and thermometer should stay in steady-state room condition for at least 30 minutes before measuring.
- Don't take a measurement while or immediately after nursing a baby.

- Don't use the non-contact thermometer in high humidity environments.
- Patients should not drink, eat, or be physically active before/while taking the measurement.
- Before placing the thermometer sensor onto the measurement area, remove hair, sweat, or dirt.
- Use an alcohol swab to carefully clean the sensor and wait for 5 minutes before taking a measurement on another patient.
- 7 short beeps will sound when the temperature is higher than 37.5 °C in order to alert the patient that he/she may have a fever.
- Always take the temperature in the same location, since the temperature readings may vary from different locations.
- In the following situations it is recommended that three temperatures in the same location be taken and the highest one taken as the reading:
 - 1) Newborn infants in the first 100 days.
 - 2) Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - 3) When you're using the thermometer for the first time until you have familiarized yourself with the instrument and obtain consistent readings.

6. Explanation of displays and operation instructions

LCD display	Operation instructions and explanation of displays	Sound and backlit
	<p>1. Measurement of body temperature</p> <p>1. Aim the temperature sensor (1) at the measurer's forehead directly (effective distance 1 to 6 cm), under shutdown state press button (2) or hold down button (2), the LCD display will display the measured forehead temperature.</p>	<p>Body temperature mode</p> <p>1. When temperature value is between 32.0°C and 37.6°C there will be one long beep and green backlit for 3 seconds.</p> <p>2. When temperature value is</p>


	<p>Holding down the button and scanning over the forehead is suggested. The system will read sequentially and show the maximum value.</p> <p>2. Measurement of milk temperature. Under shutdown state Press button (3) and release it within 4 seconds. the LCD display will display the measured milk temperature.</p> <p>3. Press button (2) or (3) under startup mode to measure again.</p> <p>Attention: There will be 7 short rapid beeps when the measured temperature exceed 37.5°C/99.5°F which is a warning for patients that he/she might have a fever.</p>	between 37.6°C/99.6°F and 42.2°C/107.9°F, there will be 7 short beeps and red backlit for 3 seconds.
	<p>The feeding bottle symbol will be visible in the left top corner when in Milk temperature mode. The °C symbol will blink when the reading is ready.</p> <p>The forehead symbol will be visible in the right top corner when in Body temperature mode. The °C symbol will blink when the reading is ready.</p>	
	<p>1. The measured value exceeds 100°C/212.0°F under milk temperature mode.</p> <p>2. The measured value exceeds 42.2°C/107.9°F under body temperature mode.</p>	3 short beeps and red backlit for 3 seconds.

	<p>1. The measured value is lower than 0°C/32.0°F under milk temperature mode. 2. The measured value is lower than 32°C/89.6°F under body temperature mode.</p>	3 short beeps and red backlit for 3 seconds.
---	---	--


Retrieval of stored readings, Storing 20 readings

LCD display	Operation instructions and explanation of displays	Sound and backlit
	Press and hold down button (3) for 4-8 seconds and LCD displays "---" with the M symbol blinking.	
 	Press button (3) again and the LCD displays the first stored reading with the M symbol blinking. Press button (3) and the LCD will display the serial number of the stored reading for 1 second and then display the measured temperature. There are 20 stored readings.	
	The LCD display will only display "---", °C/°F and the M symbol with the M symbol blinking if there are no stored readings available.	

C/F conversion

LCD display	Operation instructions and explanation of displays	Sound and backlit
 	Press and hold down button (3) for 8-12 seconds under shutdown state to enter the C/F conversion mode. Press button (3) within 5 seconds after release to change the temperature unit from °C to °F or from °F to °C	

Error message

LCD display	Operation instructions and explanation of displays	Sound and backlit
	The LCD Display will display "Er1" when the ambient temperature exceeds 40.0 °C / 104°F or drops below 10.0 °C / 50.0°F	3 short beeps and red backlit for 3 seconds.
	The LCD Display will display "ErE" if there is EEPROM data reading error or if the correcting process is not finished. Please contact your supplier.	3 short beeps and red backlit for 3 seconds.
	The LCD Display will display the low-voltage signal (non blinking) when the battery voltage is below 2.61V . Please replace the battery.	

If the thermometer is not used for 10 seconds it will automatically switch off.

7. Cleaning and storage

Use an alcohol swab or cotton swab moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the instrument in water or other cleaning liquids. After cleaning the unit should not be used until the cleaning agent is completely vaporized as the cleaning agent may affect the accuracy of measurements. Take care not to scratch the surface of the LCD display.

8. Technical specifications

Items	Standards
Model	Luvion Exact non-contact thermometer
Applicable regulations and laws	ASTM 1965 GB/T 19146-2010
Temperature units	°C / °F, adjustable
Measurement range	Body temperature mode: 32.0°C - 42.2°C / 89.6°F - 107.9°F Milk temperature mode: 0.0°C - 100.0°C / 32.0°F - 212.0°F
Precision	±0.2°C / ±0.4°F
Display resolution	±0.1°C / ±0.1°F
Latency Time	1 second
Abnormal state display	LCD Displays "L °C" if the measured temperature is below the minimum of measurement range. LCD Displays "H °C" if the measured temperature is above the maximum of measurement range.

	LCD displays "Er1" if the temperature measurement circuit is abnormal (fault of SENSOR or the temperature measurement circuit), or the ambient temperature exceeds 10°C ~ 40°C (50°F ~ 104°F). There shall be "ErC" if the calibration process is not completed or EEPROM is abnormal.
Sound and Backlit	Volume ≥50 db (the perpendicular distance from dB Volume sensor to thermometer is 10cm). Two kinds of backlit: red and green Body temperature mode: <ul style="list-style-type: none">• There will be one long beep and green backlit for 3 seconds when the temperature is between: 32.0°C/89.9°F and 37.6°C/99.6°F• There will be seven short beeps and red backlit for 3 seconds when the temperature is between 37.6°C/99.6°F and 42.2°C/107.9°F Milk temperature mode: <ul style="list-style-type: none">• There will be one long beep and green backlit for 3 seconds when the temperature is between: 0.0°C/32.0°F and 100.0°C/212.0°F
Automatic shutdown function	10s ±1s
Low-voltage display function	The product shall display low-voltage signal if the voltage is below 2.61V ±0.15V
Memory function	Memorize 20 groups of measured temperature
Current consumption	I _{stand-by} < 2 μA; I _{working} < 0.5mA; I _{Buzzer on} < 2mA; I _{Backlit} < 15mA
LED backlit specifications	Red ≥ 1.2cd/m ² Green ≥ 2.0cd/m ²
Operational conditions	ASTM 10°C - 40°C / 50°F - 104°F 15-95% RH

Type of measuring	Applicable for body temperature and milk temperature
Battery	Two 1.5V AAA batteries
Battery life	More than a 1000 times
Accuracy for clinical test	<p>The maximum allowable error for clinical test is specified in the formula below:</p> $T_{\text{error}} = \frac{(T_1 - T_{\text{ref}}) + (T_2 - T_{\text{ref}})}{2}$ <p>$\leq 0.3^\circ\text{C} / 0.6^\circ\text{F}$ (for 95%)</p> <p>T_1 and T_2 represent temperature value for the thermometer under test respectively, T_{ref} represents the constant reference temperature.</p>

It's recommended that the professional users apply technical inspection once every two years in accordance with Medical Product User Law.

9. After-sale service

The device is under warranty for one year since the date of acquisition. Application for repairing should be presented during the warranty period. The damage caused by improper use is not under warranty scope. Batteries and packaging are not under warranty scope as well.

10. Security type

 The signal indicates that the thermometer is a BF type device with internal power supply.

Please do not dispose of the product and its batteries in the household waste at the end of its useful life.

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