LUVION®

EXACT-80 NON-CONTACT

Model: HTD8813

CE 0598

HANDLEIDING 

INSTRUCTION MANUAL 

MANUAL D'INSTRUCTIONS 

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FIGURES: FIGUREN EN SYMBOLEN / FIGURES AND SYMBOLS / FIGURES ET SYMBOLES / FIGUREN UND SYMBOLE..... 4

NEDERLANDS.....6

ENGLISH14

FRANÇAIS.....26

DEUTSCH.....34

PRODUCTOVERZICHT / PRODUCT LAYOUT / DESCRIPTION DU PRODUIT

NON-CONTACT THERMOMETER


Fig. 1

DISPLAY


Fig. 2

VERKLARING VAN SYMBOLEN / EXPLANATION OF SYMBOLS / LÉGENDE DES SYMBOLES / ZEICHENERKLÄRUNG

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Hetaida Technology Co., Ltd.
Room 801 ,802, 803, 804, 901, 2# Building Scientific Research Center, Songhu Intelligent Valley, No.6 Minfu Road, Liaobu Town, Dongguan City, Guangdong Province, P.R.China

Share Info Consultant Service LLC
Repräsentanzbüro, Heerdter Lohweg 83, 40549 Düsseldorf
Tel: 0049 1767 0057022, Fax: 0049 1767 0057022
E-mail: eu-rep@share-info.cn; info@share-info.cn; will.chan@share-info.cn

**PACKAGE CONTENTS:**

- Luvion Exact-80 non-contact body thermometer
- 2 batteries 1.5V (type AAA/LR03)
- This Instruction manual

PRODUCT OVERVIEW

Luvion non-contact Exact infrared body thermometer (Fig. 1)

1. Backlight LCD
2. Battery compartment
3. ON-button (Start of measurement)
4. SET-button
5. MEMO-button
6. MODE-button
7. IR Sensor

Display (Fig. 2)

8. Surface temperature measurement mode
9. Body temperature measurement mode
10. Measurement result
11. Indicator of measurement result
12. Low battery indicator
13. Acoustic signal symbol
14. Memory display
15. Measurement unit: ° Fahrenheit
16. Measurement unit: ° Celsius
17. Memory indicator

INTRODUCTION

Thank you for purchasing the Luvion Exact-80 non-contact infrared body thermometer. Make sure you always keep your original proof of purchase. For warranty cases you will always be asked to show a copy of your proof of purchase. For warranty claims you can contact the Luvion dealer or the local Luvion representative. More information can be found at www.luvion.com

This manual contains all the information you need, to be able to use this Luvion Exact-80 Non-contact infrared thermometer safely and with pleasure. Before using the thermometer, we recommend that you read this instruction manual carefully. We ask for your extra attention for the safety instructions. Keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on this instruction manual for use. This Non-contact Infrared Body Thermometer is a medical device and can be used repeatedly and has a 5-year life expectancy.

INTENDED USE

The Luvion Exact-80 non-contact infrared body thermometer is designed to be used for intermittent measurement and monitoring of human body temperature by consumers in the home

setting and primary care setting as screening tool. Additionally the device is suitable to measure the temperature of objects in the range of 0°C to 100°C (32°F - 212°F).

SAFETY INFORMATION

- Only use the unit for the specific purpose described in the instruction manual. Any misuse will void the warranty.
- The unit is not designed for commercial use.
- Using this thermometer is no substitute for consulting a doctor.
- This product is designed to measure human body temperature on the forehead. Do not use it for any other purpose.
- Do not use the thermometer if it malfunctions or has been damaged in any matter.
- This product must be operated in a stable environment. If the ambient conditions are changeable, the sensor may mist up, causing an incorrect measurement. Clean the sensor according to the section „Care and Storage“.
- When the ambient temperature of the thermometer changes too much, such as moving the Thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C (59°F- 104°F).
- Remove batteries if equipment is not likely to be used for long time.
- This product is not waterproof, do not immerse in water or other liquid; for cleaning and disinfection, please follow the „Care and Storage“ section requirements.
- Do not touch the IR sensor with your fingers.
- Avoid taking fever measurements if measures have been taken beforehand to cool the user (e.g. cold forehead compresses, etc.)
- For the measurement of human forehead temperature, select „body“ mode; for measuring other objects, liquids, food and other temperatures, select „surface“ mode.
- Do not use near strong electrostatic or magnetic fields to avoid inaccurate measurement results.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.
- Children must be supervised to ensure that they do not play with the device.
- The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by this thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer. A control measurement using a conventional thermometer is recommended in the following cases:
 1. If the reading is surprisingly low,
 2. For new-born infants up to 100 days old,
 3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever,
 4. If you are using the thermometer for the first time or want to familiarise yourself with it.
- Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with this thermometer during or after breastfeeding, because the skin temperature may then be lower than the internal body temperature.

- The device should be kept out of the reach of children/pets.
- When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the device.
- This thermometer will not replace diagnostics in hospitals.
- Do not disassemble or modify the device. Do not let it fall down.
- Do not use this device if you think it is damaged or if you notice anything unusual. The continuous use of a damaged unit may cause injury, improper results or other serious dangers.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the „Technical Specifications“ section.
- No servicing/maintenance while the thermometer is in use.
- Never touch batteries and the patient simultaneously.
- Based on the current science and technology, no potential allergic reactions are known.
- This equipment needs to be installed and put into service in accordance with the information provided in this instruction manual.

BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before putting in the batteries!
- Remove discharged batteries from the device immediately!
- Do not use rechargeable batteries!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Insert the batteries correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the batteries from the device if it is not going to be used for an extended period!
- Keep batteries out of children's reach!
- Do not attempt to recharge the batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!

OPERATION

General notes

The ambient temperature should be stable. Keep device away from the larger flow fans, direct sunlight, air-conditioning vents and so on. Wait at least 30 seconds between each measurement. If for some reason the forehead temperature measurement cannot be realized, try to align behind the ears.

Setting the measuring mode

When the device is running, pressing the **MODE-button** (6) to switch between „body“ mode and „surface“ mode. „Body“ mode is used for measuring human body temperature, the „Surface mode“ is used to measure the surface temperature (the factory default is „body“ mode).

Perform fever measurement

1. Turn on the thermometer by pressing the **ON-button** (3). The thermometer will perform a self-test with all segments displayed for 2 seconds.
2. Immediately align the thermometer with the middle of the forehead to measure body temperature (between the eyebrows) and keep a distance of 1 - 5 cm. Alternatively you may use an alignment behind the ears. After measurement, a sound will be emitted and the temperature will be displayed. If you hear three beeps, the measurement was not successful. Press the **ON-button** (3) again and repeat the measurement.
3. Depending on the measured body temperature the backlight LCD display will change the colors:
Green means normal temperature (35.5°C- 37.3°C / 95.9°F-99.1°F),
yellow means slight fever (37.4°C-38.0°C / 99.3°F-100.4°F),
red means high fever (38.1°C-43.0°C / 100.6°F-109.4°F).

Perform surface measurement

1. Turn on the thermometer by pressing the **ON-button** (3). The thermometer will perform a self-test with all segments displayed for 2 seconds.
2. Aim the sensor at the center of the object at a distance of min. 5 to max. 15 cm. If it is covered with water, dust or dirt, remove this first to optimize the measurement result.
3. In surface mode the backlight is always green.

Recalling and erasing memory data

The last temperature taken before the thermometer powers off is automatically stored in memory (max. 50 results). With the unit switched off, press the **MEMO-button** (5) to call up the values stored in memory. If not all 50 memory locations are occupied, „---°C“ or „---°F“ may appear. If you press and hold the **MEMO-button** (5) for about 2 seconds during the memory display, „CLR“ appears in the display and the stored data is deleted. In addition, an acoustic signal sounds.

Parameter settings

The device allows some individual settings (F1 to F3): Press and hold the **SET-button** (4) when the instrument is switched on to modify these settings.

F1: Measurement unit. °C or °F.

F2: Fever alarm temperature

F3: Sound on/off

Press the **MODE-button** (6) or the **MEMO-button** (5) to adjust the value and the **SET-button** (4) to confirm and proceed to the next setting (after F3 the setting mode is completed).

Restoring the factory settings

With the device switched on, press and hold the **MODE-button** (6) until „rst“ appears on the display. 2 seconds later the factory settings for parameters F1 to F3 are restored

(F1=°C; F2=38.1°C; F3=on).



Care and Storage

- The IR sensor 7 is the most delicate part of the medisana infrared body thermometer TM A79. Protect it from dirt and damage. Clean the sensor with a cotton bud dipped in warm water or 70% isopropyl alcohol.
- Use a clean, soft cloth to clean the surface of the unit and the display. Do not use solvents or immerse the product in water or other liquids.
- Always keep the thermometer within the storage temperature range (- 20°C - +55°C or -4°F - +131°F) and humidity range (\leq 93%, non-condensing).
- It is recommended to store the thermometer in a dry, dust-free place. Do not expose the thermometer to direct sunlight, high temperature/humidity or extreme environments, as this may affect its function.

WHAT IS FEVER?

- Body temperature varies from person to person and fluctuates during the course of the day. For this reason, it is suggested to know one's normal, healthy forehead temperature to correctly determine the temperature.
- Body temperature runs approximately from 35.5°C to 37.8°C (95.9°F-100°F). To determine if one has a fever, compare the temperature detected with the person's normal temperature. A rise over the reference body temperature of 1°C (1°F) or more is generally indication of fever.
- Different measurement sites (rectal, axillary, oral, frontal, auricular) will provide different readings. Therefore it is wrong to compare the measurement taken from different sites.
- Here below are typical temperatures for adults, based on different measurement sites:
 - Rectal: 36.6°C - 38°C / 97.9°F - 100.4°F
 - Axillary: 34.7°C - 37.3°C / 94.5°F - 99.1°F

ERROR MESSAGES

Display	Problem	Solution
Body °C 	Measured body temperature outside normal limits (34.0 - 43.0°C / 93.2 - 109.4°F).	Repeat the measurement in „Body“ mode in the correct way.
Body °C 		
Body °C 	Operating temperature not correct	Establish the correct operating conditions and wait 30 minutes before operating the device
Body °C 	Distance more than 1-5 cm	The optimal distance is 1 cm
	Incorrect measuring position	Please read the instructions
	Sweat, forehead hairs or similar obstacles	Remove the obstacles
	Display flickers or turns off	Replace batteries. The device may also be damaged. Contact the service center
	Ambient temperature changes too quickly	Environmental conditions must be stable
	Battery capacity too low	Replace batteries
 (blank display)	Power is off / Batteries not inserted correctly / Batteries exhausted	Replace the batteries. Insert the batteries correctly. Turn on the thermometer

If the problems cannot be solved with these instructions or other problems occur, please contact the customer service.

DISPOSAL

This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Please remove the batteries before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Consult your municipal authority or your supplier for information about disposal.

BATTERY REPLACEMENT

If the low battery indicator (12) or nothing appears in the LCD-display (1) after a certain period of use, you must replace the batteries with two new 1.5V-batteries (type AAA / LR03).

1. Slide the lid of the battery compartment cover downwards and take the cover off.
2. Take out the batteries.
3. Put in two new 1.5V-batteries (AAA / LR03). Make sure the batteries are the right way around.
4. Put the battery compartment cover back on and move it upwards until it locks into position.

TECHNICAL SPECIFICATIONS

Display system:	Digital display (LCD)
Power supply:	3 V, 2 x AAA-batteries 1.5 V
Output range:	Body mode: 34.0°C-43.0°C (93.2°F-109.4°F)
Surface mode:	0°C-100°C (32°F-212°F)
Laboratory accuracy: (laboratory readings)	Body mode: 34.0°C- 34.9°C: ±0.3°C / 93.2°F-94.8°F: ±0.5°F 35.0°C- 42.0°C: ±0.2°C / 95.0°F-107.6°F: ±0.4°F 42.1°C- 43.0°C: ±0.3°C / 107.8°F-109.4°F: ±0.5°F
Surface mode:	±2°C (±3.6°F)
Display resolution:	0.1°C/0.1°F
Measuring time:	≤ 2 seconds
Measuring distance:	1-5 cm
Reference body site:	Axillary
Memory capacity:	Up to 50 measurement results
Automatic switch-off:	≤ 18 seconds
Battery life:	approx. 2.000 measurements
Operating temperature:	15°C to 40°C (59°F - 104°F), relative humidity: ≤85%,
Storage/ atmospheric pressure:	70-106 kPa
transport temperature:	-20°C to 55°C (-4°F - 131°F), relative humidity: ≤93%,

atmospheric pressure:	70-106 kPa
Dimensions:	approx. 138 x 95 x 40 mm
Weight:	approx. 90g without batteries

Compliance	
Item:	Compliant with:
Equipment classification	Safety Standards: EN 60601-1:2006+A1:2013, EN 60601-1-2:2015
Type of protection	Internally powered equipment (on battery power)
Degree of protection	Non Applied part
Front panel and case labeling	EN ISO15223-1:2016
Temperature	EN ISO80601-2-56:2017
Home healthcare environment	EN 60601-1-11:2015

Calculated values of the indicators according to ISO 80601-2-5

EN

Forehead mode:				
-0.042	Group A1	Group A2	Group B	Group C
Bias	0.015	0.000	-0.042	-0.040
Standard Deviation	0.123	0.121	0.121	0.137
Limits of Agreement	0.245	0.243	0.243	0.274
Clinical Repeatability (for All Group)	0.071			

Note: the above value is calculated from clinical data of HTD8818A.

Safety classification of ME EQUIPMENT	
Protection against electric shock	Internally powered ME equipment
Applied part	Non Applied part
Protection against harmful ingress of water or particulate matter	IP22
Mode of operation	Continuous operation

Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT

EMC Declaration

1. This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;
This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
2. * Caution: Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
3. *Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
4. * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacture's declaration – electromagnetic emission

The Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer of the user of the Infrared Body Thermometer should assure that it is used in such an environment.

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emission CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable

Guidance and manufacture's declaration – electromagnetic immunity

The Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of Infrared Body Thermometer should assure that it is used in such an environment.

Anti-interference detection	IEC 60601-1 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 KV Air: ± 2 , ± 4 , ± 8 , ± 15 KV	Same as the left
Electrical fast transient/burst IEC 61000-4-4	The input a.c. power ports: ± 2 KV The input d.c. power ports: ± 2 KV Signal input/output ports: ± 1 KV	Not applicable
Surge IEC 61000-4-5	Input power ports: +0.5, +1.0 KV Signal input/output: +2.0 KV	
Voltage dips IEC 61000-4-11	0.5 cycles for $> 95\%$ (sync angle (degrees): 0, 45, 90, 135, 180, 225, 270, 315) 1 cycles for $> 95\%$ UT (sync angle (degrees): 0)	
Voltage interruption IEC 61000-4-11	250 (50Hz)/300 (60Hz) cycles for $> 95\%$ UT (sync angle (degrees): 0)	
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30A/m	Same as the left
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Guidance and manufacture's declaration – electromagnetic immunity

The Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Infrared Body Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	Professional healthcare environment: 3 V/m Home healthcare environment: 10 V/m 80 MHz to 2700 MHz	Same as the left. 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Infrared Body Thermometer is used exceeds the applicable RF compliance level above, the Infrared Body Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Infrared Body Thermometer.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Guidance and manufacturer's declaration - RF wireless communication equipment immunity						
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{a)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{a)} 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{a)} 18 Hz	2	0,3	28
870						
930						
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{a)} 217 Hz	2	0,3	28
1845						
1970						
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{a)} 217 Hz	2	0,3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{a)} 217 Hz	0,2	0,3	9
5500						
5785						

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.

The current version of this instruction manual can be found under www.luvion.com

The CE marking of the device refers to the EU Directive 93/42/EEC. The device is designed for domestic use. The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the instruction manual, periodic re-adjustment is not required. If you should question the accuracy of temperature measurements, please contact the service center.

Hetaida Technology Co., Ltd.
4F,5F,6F, Comprehensive Building, Baishida High-Tech Park, Ailingkan Village,
Dalingshan Town, Dongguan City, Guangdong Province, P.R. China.

Share Info Consultant Service LLC
Repräsentanzbüro, Heerdter Lohweg 83, 40549 Düsseldorf
Tel: 0049 1767 0057022, Fax: 0049 1767 0057022
E-mail: eu-rep@share-info.cn; info@share-info.cn; will.chan@share-info.cn

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